

FitRx[™]

i-massage



Get ready to relax and wind down with **FitRx i-massage**. The combination of air compression, vibration, audio, and heat pamper and soothe. Choose from five modes—auto, comfort, dynamic, passion, and sleep for a personalized experience. Perfect for office workers, night owls, students—virtually anyone suffering from eye fatigue.

What's in the box?

- Eye Massager
- Micro-USB Charging Cable

Features:

- Wireless with rechargeable battery
- Comfortable material with adjustable strap
- Five different massage modes

Benefits:

- Soothing heat relaxes you and relieves discomfort.
- Light vibrations gently massage to ease stress and tension.
- Air compression revitalizes blood flow and mitigates pain.
- Soothing built-in music or play your own via Bluetooth.

Instructions:

1. Put the eye massager over your eyes and adjust the strap to comfortably fit around your head.
2. Press and hold the power button to turn the massager ON/OFF.
3. Press the power button to cycle through the five massage modes:

Automatic Mode: air pressure, vibration, heat

Comfort Mode: air pressure, vibration

Dynamic Mode: air pressure, vibration, heat

Passion Mode: air pressure, vibration, heat

Sleep Mode: air pressure, vibration, hot compress

4. To charge the massager, simply plug the Micro-USB end of the charging cable (included) into the massager, and the other end into a USB power source.

- A blue indicator light will flash, letting you know that the device is charging. The indicator light will turn off when it is fully charged.

NOTE: If you are wearing contacts or glasses, please remove before using this device.

NOTE: Be sure to keep your eyes closed while using the massager.

relieve. recover. relax

Playing Music:

1. While in a mode that uses music, press the power button twice to turn it ON/OFF.
2. You can play your own music through the i-massage by connecting to Bluetooth. The music can only be controlled through your device.
 - The i-massage is automatically ready to pair once you turn it on. Just open your device's Bluetooth settings and connect to "Tzumi Eye Massager".

Customer Support:

If you have any questions about using your FitRx i-massage we are here to help!
Call or email us and an experienced customer support representative will be happy to assist you:
1-855-60-TZUMI or **support@tzumi.com**.

Warnings:

- Only use while sitting or laying down. Do not use while walking or driving.
- Do not clean with water.
- Avoid using for more than 30 minutes.
- Do not disassemble.
- Do not use if device has been exposed to water. Do not operate with wet hands.
- Do not use in high temperature or humid environments, such as a sauna.
- Avoid exposure to direct sunlight while using this device.
- Massager could get hot, be careful if you have sensitive eyes.
- Consult your doctor if you have sensitive skin.

Disclaimer:

Read the instructions carefully before using the eye massager. The use of this massager is subject to following its instructions closely. tzumi is not liable for any damage or loss caused by inappropriate use of the device.

If you are: pregnant, diabetic, have a pacemaker, recovering from surgery, suffer from migraines and/or neuropathic disorders, disc protrusion and/or spinal inflammation, do not use this massager without the approval of your doctor. If you feel pain or discomfort while using the massager, stop using it immediately and consult your doctor.

**For products and innovations visit
www.FitRxRecovery.com**



FitRx i-massage. Distributed by tzumi Inc. NY, NY, 10016. ©Copyright 2021 tzumi Inc. All Rights Reserved. Made In China. All other trademarks and trade names are the property of their respectful owners. All pictures shown are for illustration purpose only. Actual product may vary due to product enhancement. Subject to change without prior notice.